



# Recombinant Moloney murine sarcoma virus Gag polyprotein (gag)

<b>Product Code</b>	CSB-MP341338MHO
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P32594
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Moloney murine sarcoma virus (strain ts110) (MoMSV)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	GQVTITPLS LTLDHWKDVE RIAHNQSVDV KKRRWVTFCS AEWPTFNVGW PRDGTFNRLD ITQVKIKVFS PGPHGHPDQV PYIVTWEALA FDPWPVKPF VHPKPPPLL PSAPSLPLEP PLSTPPQSSL YPALTPSLGA KPKPQVLSDS GGPLIDLLTE DPPPYPDPRP PPSDRDGDGSG EATPAGEAPD PSPMASRLRG RREPPVADST TSQAFPLRTG GNGQLQYWPF SSSDLYNWKS NNPSFSEDPG KLTALIESVL ITHQPTWDDC QQLLGTLLTG EEKQRVLEA RKAVRGDDGR PTQLPNEVDA AFPLERPDWE YTTQAGR NHL VHYRQLLIAG LQNAGRSPTN LAKVKGITQG PNESPSAFLE RLKEAYRRYT PYDPEDPGQE TNVSMSFIWQ SAPDIGR KLE RLEDLRN KTL GDLVREAERI FNKRETPEER EERIRREREE KEERHAPKLP WLFLIISP
<b>Source</b>	Mammalian cell
<b>Target Names</b>	gag
<b>Protein Names</b>	Recommended name: Gag polyprotein Alternative name(s): Core polyprotein Cleaved into the following 3 chains: 1. Matrix protein p15 Short name= 2. MA 3. RNA-binding phosphoprotein p12 Alternative name(s): pp12 Capsid protein
<b>Expression Region</b>	2-468
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.