



Recombinant Escherichia coli Uncharacterized protein yjcS (yjcS)

Product Code	CSB-YP341349ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P32717
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MNNSRLFRLS RIVALTAAS GMMVNTANAK EEAKAATQYT QQVNQNYAKS LPFSDRQDFD DAQRGFIAPL LDEGILRDAN GKVYYRADDY KFDINAAPE TVNPSLWRQS QINGISGLFK VTDKMYQVRG QDISNITFVE GEKGIIVIDP LVTPPAKAA LDLYFQHRPQ KPIVAVIYTH SHTDHYGGVK GIISEADVKS GKVQVIAPAG FMDEAISENV LAGNIMSRA LYSYGLLLPH NAQGNVGNGL GVTLATGDPS IIAPTKTIVR TGEKMIIDGL EFDLMTPGS EPAEMHFYI PALKALCTAE NATHTLHNFY TLRGAKTRDT SKWTEYLNEL LDMWGNDAEV LFMPHTWPVW GNKHINDYIG KYRDTIKYIH DQTLHLANQG YTMNEIGDMI KLPPALANNW ASRGYYGSVS HNARAVYNFY LGYYDGNPAN LHPYGQVEMG KRYVQALGGS ARVINLAQEA NKQGDYRWSA ELLKQVIAAN PGDQVAKNLQ ANNFEQLGYQ AESATWRGFY LTGAKELREG VHKFSHGTTG SPDTIRGMSV EMLFDFMAVR LDSAKAAGKN ISLNFNMSNG DNLNLTNDS VLNYRKTLPQ QADASFYISR EDLHAVLTGQ AKMADLVKAK KAKIIGNGAK LEEIIACLDN FDLWVNIVTP N
Source	Yeast
Target Names	yjcS
Protein Names	Recommended name: Uncharacterized protein yjcS
Expression Region	1-661
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.