



Recombinant Arabidopsis thaliana Serine carboxypeptidase-like 49 (SCPL49)

Product Code	CSB-YP341361DOA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P32826
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	FPNILADG GPTVDDLGHG AGYYKLPKSR GASMFYFFFE SRNKKDAPVV IWLTTGGPGCS SELAVFYENG PFKITSNMSL AWNEYGWDQV SNLLYVDQPV GTGFSYTTDK SDIRHDETGV SNDLYDFLQA FFAEHPKLAK NDFYITGESY AGHYIPAFAS RVHKGKANE GVHINLKGFA IGGLTDPAL QYPAYPDYAL EMGLITQKEH DRLEKIVPLC ELSIKLCGTD GTTSCCLASYL VCNSLFSGVM SHAGGVNYYD IRKKCVGSLC YDFSNMEKFL NLQSVRKS LG VGDIDFVSCS TSVYQAMLVD WMRNLEVGIP TLLEDGISLL VYAGEYDLIC NWLGNRWWN AMEWSGKTNF GAAKEVPFIV DGKEAGLLKT YEQLSFLKVR DAGHMVPMDQ PKAALKMLKR WMENSLIEDA TVTVAAQGGG ELVAQM
Source	Yeast
Target Names	SCPL49
Protein Names	Recommended name: Serine carboxypeptidase-like 49 EC= 3.4.16.-
Expression Region	83-516
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.