



Recombinant Human Tripartite motif-containing protein 38 (TRIM38)

Product Code	CSB-BP024475HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O00635
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	<p>MASTTSTKKM MEEATCSICL SLMTNPVSIN CGHSYCHLCI TDDFKNPSQK QLRQETFCCP QCRAPFHMDL LRPNKQLGSL IEALKETDQE MSCEEHGEQF HLFCEDEGQL ICWRCERAPQ HKGHTTALVE DVCQGYKEKL QKAVTKLKQL EDRCTEQKLS TAMRITKWKE KVQIQRQKIR SDFKNLQCFL HEEKSYLWR LEKEEQQTLR RLRDYEAGLG LKSNELKSHI LELEEKQCGS AQKLLQNVND TLRSWAVKL ETSEAVSLEL HTMCNVSKLY FDVKKMLRSH QVSVTLDPDT AHHELILSED RRQVTRGYTQ ENQDTSSRRF TAFPCVLGCE GFTSGRRYFE VDVGEETGWD LGVCMENVQR GTGMKQEPQS GFWTLRLCKK KGYVALTSPV TSLHLHEQPL LVGIFLDYEA GVVSFYNGNT GCHIFTFPKA SFSDTLRPYF QVYQYSPLFL PPPGD</p>
Source	Baculovirus
Target Names	TRIM38
Protein Names	Recommended name: Tripartite motif-containing protein 38 Alternative name(s): RING finger protein 15 Zinc finger protein RoRet
Expression Region	1-465
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Target Details	This protein is a member of the tripartite motif (TRIM) family. The TRIM motif includes three zinc-binding domains, a RING, a B-box type 1 and a B-box type 2, and a coiled-coil region. The function of this protein has not been identified.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.