



Recombinant Cat Granulocyte colony-stimulating factor (CSF3)

Product Code	CSB-MP006049CA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O02708
Product Type	Recombinant Protein
Immunogen Species	Felis catus (Cat) (Felis silvestris catus)
Purity	≥85% (SDS-PAGE)
Sequence	TPLGPTSSLP QSFLKCLEQ VRKVQADGTA LQERLCAAHK LCHPEELVLL GHALGIPQAP LSSCSSQALQ LTGCLRQLHS GLFLYQGLLQ ALAGISPELA PTLDMLQLDI TDFAINIWQQ MEDVGMAPAV PPTQGTMPPTF TSAFQRRAGG TLVASNLQSF LEVAYRALRH FTKP
Source	Mammalian cell
Target Names	CSF3
Protein Names	Recommended name: Granulocyte colony-stimulating factor Short name= G-CSF
Expression Region	21-194
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This protein is a cytokine that controls the production, differentiation, and function of granulocytes. The active protein is found extracellularly. Three transcript variants encoding three different isoforms have been found for this gene.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.