



Recombinant Uncharacterized calcium-binding protein C07A9.5 (C07A9.5)

Product Code	CSB-EP341515CXY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P34316
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	MIEKKIIADK FNSNFIWIFV DAADKYPKQM NIKANQEHII YGEKCRLVGF LETESPSTGR FSTNHSTWTH SRPTVSGERF THSFKTLASR MERLDQTISD RIRIVDGEIS RHDNEKICE QAACTLKTCN NVIVDVKKEL EMLHSLRLED QKSLMNLIA KVQKSGMIQD LEKWRSLMES RYIFDNTFSS ATPHGVLVEM CQCLELMSSM LRSVEQSIAD RNHNGVTEKQ LHEFKLAFDY FDQEKNGWLD YEHFELCLKS QGYNFSVEST LKETMSLLDP STYEKHDYMR YMVKHETTNI LDDHSAIEDA LKNLDARKIS DSMSRKEAEF FMSKIAKKA E TSSDQVCLEY KDFVNSFY
Source	E.coli
Target Names	C07A9.5
Protein Names	Recommended name: Uncharacterized calcium-binding protein C07A9.5
Expression Region	1-358
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.