



Recombinant Bovine CCAAT/enhancer-binding protein beta (CEBPB)

Product Code	CSB-BP005181BO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O02755
Product Type	Recombinant Protein
Immunogen Species	Bos taurus (Bovine)
Purity	>85% (SDS-PAGE)
Sequence	<p> MQRLVVWDPV CLPLPPPPPA FKSMEVANFY YEADCLAAAY GGKAAPAAPP ADRPGPRPPT GELGSIGEHE RAIDFSPYLE PLGAPQAPAP TTASDTFEAA PSAPAPVPAS SGQHHDFLSD LFSDDYGGKN CKKAAEYGYV SLGRLGAAKG ALHPGCFAPL HPPPPPPPPP AELKAEPGFE PADCKRKEEA GAPGGGAAGM AAGFPYALRA YLG YQAVPSG SSGSLSTSSS SSPPGTPSPA DAKATPAAAA CYAGAAPAPS QVKSKAKKTV DKHSDEYKIR RERNNI AVRK SRDKAKMRNL ETQHKVLELT GENERLQKKV EQLSREVSTL RNLFKTLPEP LLASSGHC </p>
Source	Baculovirus
Target Names	CEBPB
Protein Names	Recommended name: CCAAT/enhancer-binding protein beta Short name=C/EBP beta
Expression Region	1-348
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Target Details	<p>The protein encoded by this intronless gene is a bZIP transcription factor which can bind as a homodimer to certain DNA regulatory regions. It can also form heterodimers with the related proteins CEBP-alpha, CEBP-delta, and CEBP-gamma. The encoded protein is important in the regulation of genes involved in immune and inflammatory responses and has been shown to bind to the IL-1 response element in the IL-6 gene, as well as to regulatory regions of several acute-phase and cytokine genes. In addition, the encoded protein can bind the promoter and upstream element and stimulate the expression of the collagen type I gene.</p>
Reconstitution	<p>We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.</p>



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.