



# Recombinant Bovine CCAAT/enhancer-binding protein beta (CEBPB)

<b>Product Code</b>	CSB-EP005181BO-B
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	O02755
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Bos taurus (Bovine)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	<p>                     MQRLVVWDPV CLPLPPPPPA FKSMEVANFY YEADCLAAAY GGKAAPAAPP                      ADRPGPRPPT GELGSIGEHE RAIDFSPYLE PLGAPQAPAP TTASDTFEAA                      PSAPAPVPAS SGQHHDFLSD LFSDDYGGKN CKKAAEYGYV SLGRLGAAKG                      ALHPGCFAPL HPPPPPPPPP AELKAEPGFE PADCKRKEEA GAPGGGAAGM                      AAGFPYALRA YLGYQAVPSG SSGSLSTSSS SSPPGTPSPA DAKATPAAAA                      CYAGAAPAPS QVKSKAKKTV DKHSDEYKIR RERNNIAVRK SRDKAKMRNL                      ETQHKVLELT GENERLQKKV EQLSREVSTL RNLFKTLPEP LLASSGHC                 </p>
<b>Source</b>	E.coli
<b>Target Names</b>	CEBPB
<b>Protein Names</b>	Recommended name: CCAAT/enhancer-binding protein beta Short name=C/EBP beta
<b>Expression Region</b>	1-348
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Target Details</b>	<p>The protein encoded by this intronless gene is a bZIP transcription factor which can bind as a homodimer to certain DNA regulatory regions. It can also form heterodimers with the related proteins CEBP-alpha, CEBP-delta, and CEBP-gamma. The encoded protein is important in the regulation of genes involved in immune and inflammatory responses and has been shown to bind to the IL-1 response element in the IL-6 gene, as well as to regulatory regions of several acute-phase and cytokine genes. In addition, the encoded protein can bind the promoter and upstream element and stimulate the expression of the collagen type I gene.</p>
<b>Reconstitution</b>	<p>We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.</p>



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.