



Recombinant Autophagy-related protein 13 homolog (atg-13)

Product Code	CSB-MP341524CXY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P34379
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MVNEYDTYNK WLKFFSVRMV QSIIQSRLGD EIESKCV PYS ENAVDWFNMR IDELGEISAY LKSNISYPP VGTLTLEFL YTPSGQLLPL EAWILSSSEG TDDCSRNELY HDMSTLLRSA IVSARMTPMH RLYVKKQHLE TFVIMYRVFE NDISSDMGKG KKTRKIGELV SKFGNISLDL HYRTSMHFEE PEIAPVTPVE DVEEIIDGGE KTVVDENIQT RTVSECVPPIA DAKKRKASGS VESATSAGSS TSREAAPRFI LGQSTSSSEDS RHSDVQNSYE EDHKPSLADL RNHSFPFVNL LQSAYNPANG TTKNSSSTCL NSPKSTPEDK EPTIEKVAES FRAAKIDEVV FEEDEDEELP LDSMELSEDS FVHFNQLSDF GGAPSLGNEL GDYLKQLKTA PDMTESGDID ICNMDLKTEL EKISSQTANF NNFLKHVNSF SDE
Source	Mammalian cell
Target Names	atg-13
Protein Names	Recommended name: Autophagy-related protein 13 homolog
Expression Region	1-443
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.