



Recombinant Uncharacterized protein F02A9.4b (F02A9.4)

Product Code	CSB-EP341525CXY-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P34384
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MFRHVAQNLG SRNTSIQSYR LLRTRWERGY LKDLYHRRQI LGADPAISRS SYPNWDYDSE LYAFGHRIGA PEISEDQYIK ALTNESFFQR ADVEENVESA QPGAIEVGTET NAELVKRGEQ NLSIWLKRYL RFHLAKAPEE LIEAIDSHLL DDECLAGIAS HLGIDHLVRT KEFPISQESS ADAFRALAGV FSDEKVKNLV IDFIVPQLVD IDFADIYPLA DPLAVVTDLL KSNQVTEIEP RVLRSAGENS AEPIYVVAIY ADKLNKVGQS AGESLAIAMD MAAREALLRL WDTSEKVLV FGDRAANVPL ERYSQPNFSL SQKCSPGTNT NIIDHSQPVE SPDNLIEAVL RYRNVVDAEV GKSQYTKRLRH KFSRGLAKR SFRYLVKPKP YTVA
Source	E.coli
Target Names	F02A9.4
Protein Names	Recommended name: Uncharacterized protein F02A9.4b
Expression Region	1-394
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.