



Recombinant Uncharacterized protein K01B6.3 (K01B6.3)

Product Code	CSB-EP341538CXY-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P34491
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	MDDDFEEVGG RFRRFRRLK KDWSQLWSDL SDPKQMRFPW PINTPLEVYE QVGAERFDPS SHRYQMPFFW KTLAVPGTPS WALLVCIIES TFALICFLAV MLHYAIFLPW CEVYSKPLMV FFLTALQYAI FYSFKVIFMI SIVERNARLL RLQLFFQYAT CVLLLLDAAF ALAADFGGYN EELIYCNKNP LLIRFVAIIS LIFLVMQFLL RIITVQVYNF MWDVRKFRKS LNNSKWRYRK RVHFTYCSIM QEDFKNERLQ NKTKSVESRF RDGQEEILKQ IQRKQNVNTNI SIAPDEPFDE AQLPQNARSI STLSTTPSNS ATVSSPLGKR KHHTVSVSSK KMRDNEGKKV KVSRRKPGGI KVQLEVDYET ARLLFSPKKN GRIPNVSVEE VDMEFDEDDD SHPEEHEPLV
Source	E.coli
Target Names	K01B6.3
Protein Names	Recommended name: Uncharacterized protein K01B6.3
Expression Region	1-410
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.