



# Recombinant Uncharacterized protein K01B6.3 (K01B6.3)

<b>Product Code</b>	CSB-MP341538CXY
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P34491
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Caenorhabditis elegans
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MDDDFEEVGG RFRRFRRLK KDWSQLWSDL SDPKQMRFPW PINTPLEVYE QVGAERFDPS SHRYQMPFFW KTLAVPGTPS WALLVCIIES TFALICFLAV MLHYAIFLPW CEVYSKPLMV FFLTALQYAI FYSFKVIFMI SIVERNARLL RLQLFFQYAT CVLLLLDAAF ALAADFGGYN EELIYCNKNP LLIRFVAIIS LIFLVMQFLL RIITVQVYNF MWDVRKFRKS LNNSKWRYRK RVHFTYCSIM QEDFKNERLQ NKTKSVESRF RDGQEEILKQ IQRKQNVNTNI SIAPDEPFDE AQLPQNARSI STLSTTPSNS ATVSSPLGKR KHHTVSVSSK KMRDNEGKKV KVSRRKPGGI KVQLEVDYET ARLLFSPKKN GRIPNVSVEE VDMEFDEDDD SHPEEHEPLV
<b>Source</b>	Mammalian cell
<b>Target Names</b>	K01B6.3
<b>Protein Names</b>	Recommended name: Uncharacterized protein K01B6.3
<b>Expression Region</b>	1-410
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.