



Recombinant Uncharacterized protein T16H12.1 (T16H12.1)

Product Code	CSB-EP341548CXY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P34564
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MSLRLTSSLN QNKTTNVNGA ISRCLWRTET RCSTGYGSKS GFAEEVGRKK LEIGRLELKW ERKLDSDLF VRRKFEKILF PNGNKSSRKK ANLDKNRRIA TKVANRFELS LNQPVHFSQT NPKSTPEPPC TSSSGAGDCH ENLPADGYPI ENDDGMGEPA RMLGKIKKPR SRLECIFCKD SHTRSIVNCP NIRKISDRTN ILIVQGRKID NTEKQRKAKK TDVRVGRKFD ALRNFIERNI EEVARDEPIS EQKVNHRSEF SEQIIFGDYE RGPTFLTGNL VSSGLSTFTM DSEVSNDSFG IPGTMPSPKY LVSPISNPTF SNKCTRLIRQ CHQLHSRNCA YFVEMSIFQ
Source	E.coli
Target Names	T16H12.1
Protein Names	Recommended name: Uncharacterized protein T16H12.1
Expression Region	1-349
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.