



# Recombinant Uncharacterized protein T16H12.1 (T16H12.1)

<b>Product Code</b>	CSB-MP341548CXY
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P34564
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Caenorhabditis elegans
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSLRLTSSLN QNKTTNVNGA ISRCLWRTET RCSTGYGSKS GFAEEVGRKK LEIGRLELKW ERKLDSDLF VRRKFEKILF PNGNKSSRKK ANLDKNRRIA TKVANRFELS LNQPVHFSQT NPKSTPEPPC TSSSGAGDCH ENLPADGYPI ENDDGMGEPA RMLGKIKKPR SRLECIFCKD SHTRSIVNCP NIRKISDRTN ILIVQGRKID NTEKQRKAKK TDVRVGRKFD ALRNFIERNI EEVARDEPIS EQKVNHRSEF SEQIIFGDYE RGPTFLTGNL VSSGLSTFTM DSEVSNDSFG IPGTMPSPKY LVSPISNPTF SNKCTRLIRQ CHQLHSRNCA YFVEMSIFQ
<b>Source</b>	Mammalian cell
<b>Target Names</b>	T16H12.1
<b>Protein Names</b>	Recommended name: Uncharacterized protein T16H12.1
<b>Expression Region</b>	1-349
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.