



Recombinant Uncharacterized protein T20G5.4 (T20G5.4)

Product Code	CSB-MP341550CXY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P34577
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	MTVSYNQSVA TSRPWTFLAL IFRWRGSVWS AIWIQYSVWL GLYFLVSAIY RFILSAYQQQ IFVRLVDYVN SRMSYVPLDW MLGFFIAGVL RRFWYLYDII GFIDNIACST ATYIRGDSER AKQYRRNIIR YCELTQVLIF RDLSMKARKR FPTLDTVAAA GFMMPEKAN FDLIQYNYNK YFLPFNWAWA LVYNARKEGL IEGDYYVTVI SEDIKKFRTG LAWVCNYDWV PLPIIYPTIV CLAVHMYFFV GILARQYVKG SEIDPDMIDL VFPFMTSIQF VFYMGWLKVG ELLNPWGED PDDFETNMLI DRNLAMGLKI VDEGYDKTPR LEKDAFWDDT WVPLYSEASA HEKRYHQKQG SLAHIKIGRS VSQVRMVPRD GRRASIVKER IVNVKPGSGI GDILRPSSLL NLMKHASSSR SLERQRSPGS FRMETLTPGS PTNTPIEPID KIDKKK
Source	Mammalian cell
Target Names	best-18
Protein Names	Recommended name: Uncharacterized protein T20G5.4
Expression Region	1-456
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.