



Recombinant Tubulin alpha-2 chain (tba-2)

Product Code	CSB-BP341566CXY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P34690
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	MREVISIHVG QAGVQIGNAC WELYCLEHGI QPDGTMPTQS TNEGESFTTF FSDTGSGRYV PRSIFVDLEP TVVDEIRTGT YKKLFHPEQM ITGKEDAANN YARGHYTVGK ELIDTVLDRI RRLADNCSGL QGFFVFHSFG GGTGSGFTSL LMERLSVDYG KSKLEFSIY PAPQVSTAVV EPYNSILTTH TTLEHSDCAF MVDNEAIYDI CRRNLDVERP SYTNLNRIIS QVVSSITASL RFDGALNVDL NEFQTNLVPY PRIHFPLAAY TPLISAEKAY HEALSVDIT NSCFEPANQM VKCDPRHGKY MAVCLLYRGD VVPKDVNTAI AAIKTKRTIQ FVDWCPTGFK VGINYQPPTV VPGGDLAKVP RAVCMLSNTT AIAEAWSRLD YKFDLMYAKR AFVHWYVGEG MEEGEFTEAR EDLAALEKDY EEVGADSNEG GEEEGEEY
Source	Baculovirus
Target Names	tba-2
Protein Names	Recommended name: Tubulin alpha-2 chain
Expression Region	1-448
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.