



Recombinant *Synechococcus elongatus*

Uncharacterized protein CcmO (ccmO)

Product Code	CSB-YP341598FPY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P46205
Product Type	Recombinant Protein
Immunogen Species	<i>Synechococcus elongatus</i> (strain PCC 7942) (<i>Anacystis nidulans</i> R2)
Purity	≥85% (SDS-PAGE)
Sequence	G TADMMLKSAD VTLIGYEKTG SGFCTAIIRG GYADIKLALE AGVATARQFE QYVSSTILPR PQGNLEAVLP ISRRLSQEAM ATRSHQNVGA IGLIETNGFP ALVGAADAML KSANVKLIC Y EKTGSGGLCTA IVQGTVSNVT VAVEAGMYAA ERIGQLNAIM VIPRPLDDL DLSLPEPQSDS EAAQPLQLPL RVREKQPLLE LPELERQPIA IEAPRLAEE RQSALELAQE TPLAEPLLEP NPRDDQ
Source	Yeast
Target Names	ccmO
Protein Names	Recommended name: Uncharacterized protein CcmO
Expression Region	30-276
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.