



Recombinant Triticum aestivum Mitochondrial outer membrane porin (VDAC1)

Product Code	CSB-BP341605TQN
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P46274
Product Type	Recombinant Protein
Immunogen Species	Triticum aestivum (Wheat)
Purity	>85% (SDS-PAGE)
Sequence	MGGPGLYSGI GKKAKDLLYR DYQTDHKFTL TTYTANGPAI TATSTKKADL TVGEIQSQIK NKNITVDVKA NSASNVITTI TADDLAAPGL KTILSFAVPD QKSGKVELQY LHDYAGINAS IGLTANPVVN LSGAFGTSAL AVGADVSLDT ATKNFAKYNA ALSYTNQDLI ASLNLNKGKGD SLTASYHHIV EKSGTAVGAE LTHSFSSNEN SLTFGTQHTL DPLTLVKARI NNSGKASALI QHEFMPKSLC TISAEVDTKA IEKSSKVGIA IALKP
Source	Baculovirus
Target Names	VDAC1
Protein Names	Recommended name: Mitochondrial outer membrane porin Alternative name(s): Voltage-dependent anion-selective channel protein Short name= VDAC
Expression Region	1-275
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.