



Recombinant Probable leucine carboxyl methyltransferase 1 (B0285.4)

Product Code	CSB-EP341627CXY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P46554
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MDSEAVSSDS HVAAAIATRR RSNSVSDDYS VQRTNDDATQ CKYFATQKGY WKDEFISRFA NSSSNVSEAR RFPEISMGYW ARTAAIEKYV RDFLNEFDGN AQVVS LGCGF DTLFWRLVSS GAKLVKYVEV DFSSVTSKKI RHILKPIGPN SVDLKKSFES DAVVSHHADL HAGNYHLIGA DLRQANELDQ KLATCQLSHD IPTIFIAECV LVYMSADSST ALLKQIVSQF KQPAFVNYEQ FRTSDAFTKV MEQNLGDRGI QLHGLEMCES AEKQEERFRN AGFKEVKVMD MNQIFNNFLD QKEVSRIREI EMLDEMELLQ QLFAHYCVVS ARI
Source	E.coli
Target Names	B0285.4
Protein Names	Recommended name: Probable leucine carboxyl methyltransferase 1 EC=2.1.1.233 Alternative name(s): [Phosphatase 2A protein]-leucine-carboxy methyltransferase 1
Expression Region	1-333
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.