



Recombinant *Saccharomyces cerevisiae* Cop9 signalosome complex subunit 12 (CSN12)

Product Code	CSB-MP341670SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P47130
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MDVDIGCYFE EKRYDDKLLD FIRYDVKTPK KTKYILQRPT ATDEESVRLQ RFYQLGVDLK LKYSKRRSLK KQGRIKNATE ELLRLANEQL KLFNRIVERE TNWIIYPLWV MAKQLIRLAN ESSELNKDSI EECGRTIHRF FTICLNDRNP RLNENKKIGC YMFANLEFSI YHRLSNKDMI KNLVKVLESR VNARDIPPLN KSLAMEHKSQ VVLYNYYLGQ YYGCLENDHE RGGFFHLNEAL LQCPMLYVES TGKFLVQGQM EKIMILLVPL ALLTKRLYPH WDHPVIAGVI TRSKRLSQVY PTLVRSVISG NLSLYEATAA SHERFFLSQG LHVVITLLRE VVFTRLVQRC WQWGNDRKSI MPLKILLATK QHDSSANEDE EEQLDALECR LASAIASGLL RAYLSHSNRC IVFSKKEFPF HSK
Source	Mammalian cell
Target Names	CSN12
Protein Names	Recommended name: Cop9 signalosome complex subunit 12
Expression Region	1-423
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.