



Recombinant Chicken Heat shock protein beta-1 (HSPB1)

Product Code	CSB-EP010833CH-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q00649
Product Type	Recombinant Protein
Immunogen Species	Gallus gallus (Chicken)
Purity	>85% (SDS-PAGE)
Sequence	MAERRVPFTF LTSPSWEPFR DWYHGSRLFD QSFGMPHIPE DWYKWPSGSA WPGYFRLPS ESALLPAPGS PYGRALSELS SGISEIRQSA DSWKVTLDVN HFAPEELVVK TKDNIVEITG KHEEKQDEHG FISRCFTRKY TLPPGVEATA VRSSLSPDGM LTVEAPLPKP AIQSSEITIP VTVEAKKEEP AKK
Source	E.coli
Target Names	HSPB1
Protein Names	Recommended name: Heat shock protein beta-1 Short name= HspB1 Alternative name(s): 25 kDa IAP Actin polymerization inhibitor Heat shock 25 kDa protein Short name= HSP 25 Heat shock 27 kDa protein Short name= HSP 27
Expression Region	1-193
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Target Details	This protein is induced by environmental stress and developmental changes. The encoded protein is involved in stress resistance and actin organization and translocates from the cytoplasm to the nucleus upon stress induction. Defects in this gene are a cause of Charcot-Marie-Tooth disease type 2F (CMT2F) and distal hereditary motor neuropathy (dHMN).
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.