



# Recombinant Mouse Insulin-like 3 (Insl3), partial

<b>Product Code</b>	CSB-BP011747MO
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	O09107
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Source</b>	Baculovirus
<b>Target Names</b>	Insl3
<b>Protein Names</b>	Recommended name: Insulin-like 3 Alternative name(s): Leydig insulin-like peptide Short name= Ley-I-L Relaxin-like factor Cleaved into the following 2 chains: 1. Insulin-like 3 B chain 2. Insulin-like 3 A chain
<b>Expression Region</b>	-
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Partial
<b>Target Details</b>	This protein is an insulin like hormone produced mainly in gonadal tissues in males and females. Studies of the mouse counterpart suggest that this gene may be involved in the development of urogenital tract and female fertility. It may also act as a hormone to regulate growth and differentiation of gubernaculum, and thus mediating intra-abdominal testicular descent. The mutations in this gene may lead to, but not a frequent cause of, cryptorchidism.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.