



Recombinant *Saccharomyces cerevisiae* Adenosine kinase (ADO1)

Product Code	CSB-BP342064SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P47143
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MTAPLVVLGN PLLDFQADVT AEYLAKYSLK ENDAILVDAK SGDAKMAIFD ELLQMPETKL VAGGAAQNTA RGAAYVLGAG QVVYFGSVGK DKFSERLLNE NEKAGVKSMY QVQNDIGTGK CAALITGHNR SLVTDLGAAN FFTPDHLDKH WDLVEAAKLF YIGGFHLTVS PDAIVKLGQH AKENSKPFVL NFSAPFIPHV FKDALARVLP YATVIIANES EAEAFCDAFQ LDCANTDLEA IAQRIVKDSP VEKTVIFTHG VEPTVVVSSK GTSTYPVKPL DSSKIVDTNG AGDAFAGGFM AGLTKGEDLE TSIDMGQWLA ALSIQEVGPS YPSEKISYSK
Source	Baculovirus
Target Names	ADO1
Protein Names	Recommended name: Adenosine kinase EC= 2.7.1.20
Expression Region	1-340
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.