



Recombinant *Drosophila melanogaster* Protein atonal (ato)

Product Code	CSB-BP342195DLU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P48987
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	≥85% (SDS-PAGE)
Sequence	MSSSEIYRYYY YKTS EDLQGF K TAAAE P YFN PMAAYNPGVT HYQFNGNTLA SSSNYLSANG FISFEQASSD G WISSSPASH RSESPEYVDL NTMYNGGCNN MAQNQQYGM I MEQSVVSTAP AIPVASPPAV EVMGSSNVGT CKTIPASAGP KPKRSYTKKN QPSTTATSTP TAAAE SSV NLYTEEFQNF DFDNSALFDD SVEDDEDLML FSGGEDFDGN DGSFDLADGE NQDAAAGGSG KKRRGKQITP VVKRKRRLAA NARERRRMQN LNQAFDRLRQ YLPCLGNDRQ LSKHETLQMA QTYISALGDL LR
Source	Baculovirus
Target Names	ato
Protein Names	Recommended name: Protein atonal
Expression Region	1-312
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.