



Recombinant Arabidopsis thaliana Glutamyl-tRNA reductase 2, chloroplastic (HEMA2)

Product Code	CSB-BP342213DOA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P49294
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	ASISAL EQLKTSIDR YTKERSSIVV IGLSIHTAPV EMREKLAIPE AEWPRAIAEL CGLNHIEEAA VLSTCNRMEI YVLALSQHRG VKEVTEWMSK TSGIPVSEIC QHRFLLYNKD VTQHIFEVSA GLDSLVLGEG QILAQVKQVV KVGQGVNGFG RNISGLFKHA ITVGKRV RTE T NIAAGAVSV SSAAVELALM KLPESSHASS ARMLVVGAGK MGKLVIKHLV AKGCTKMVVV NRSEEKVA AV RNEMPPGVEI IYKPLDEMLS CAAEADVFT STASETPLFL KEQVETLPPV RDARLFVDIS VPRNVGSCVA EIDGTRVFNV DDLKEVVAAN KEDRVRKAMD AQAIITDESK HFEAWRDSLE TVPTIKKLRG YTERIIAAEI EKSLPKMGID MNKKMRKTVD DLIRGIVNKL LHGPMQHLRC DGND SRTLSE TLDNMQALNR MYGLDAEILE EKIRAKVEKK
Source	Baculovirus
Target Names	HEMA2
Protein Names	Recommended name: Glutamyl-tRNA reductase 2, chloroplastic Short name= GluTR EC= 1.2.1.70
Expression Region	65-530
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.