



Recombinant Arabidopsis thaliana Dihydroflavonol-4-reductase (DFRA)

Product Code	CSB-EP342342DOA-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P51102
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MVSQKETVCV TGASGFIGSW LVMRLLEGRY FVRATVRDPG NLKKVQHLLD LPNAKLLTL WKADLSEEGS YDDAINGCDG VFHVATPMDF ESKDPENEVI KPTVNGMLGI MKACVKAKTV RRFVFTSSAG TVNVEEHQKN VYDENDWSDL EFIMSKKMTG WMYFVSKTLA EKAAWDFAEK KGLDFISIIP TLVVGPFITT SMPPSLITAL SPITRNEAHY SIIRQGQYVH LDDLCNAHIF LYEQAAAKGR YICSSHDATI LTISKFLRPK YPEYNVPSTF EGVNENLKSI EFSSKKLTDM GFNFKYSLEE MFIESIETCR QKGFLPVSL YQSISEIKTK NENIDVKTGD GLTDGMKPCN KTETGITGER TDAPMLAQQM CA
Source	E.coli
Target Names	DFRA
Protein Names	Recommended name: Dihydroflavonol-4-reductase Short name= DFR EC= 1.1.1.219 Alternative name(s): Dihydrokaempferol 4-reductase Protein TRANSPARENT TESTA 3
Expression Region	1-382
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.