



# Recombinant Human Mitochondrial ribonuclease P protein 3 (KIAA0391)

<b>Product Code</b>	CSB-BP012187HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	O15091
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	FSLKT MSPQNTKATN LIAKARYLRK DEGSNKQVYS VPHFFLAGAA KERSQMNSQT EDHALAPVRN TIQLPTQPLN SEEDKDKED LKENTGKTSF ESWIISQMAG CHSSIDVAKS LLAWVAAKNN GIVSYDLLVK YLYLCVFHMQ TSEVIDVFEI MKARYKTLEP RGYSLIRGL IHSDRWREAL LLEDIKKVI TPSKKNYNDQ IQGALLHQDV NTAWNLYQEL LGHDIVPMLE TLKAFFDFGK DIKDDNYSNK LLDILSYLRN NQLYPGESFA HSIKTFWESV PGKQWKGQFT TVRKSGQCSG CGKTIESIQL SPEEYECLKG KIMRDVIDGG DQYRKTTTPQE LKRFFENFIKS RPPFDVVIDG LNVAKMFPKV RESQLLLNVV SQLAKRNLRL LVLGRKHMLR RSSQWSRDEM EEVQKQASCF FADDISEDDP FLLYATLHSG NHCRFITRDL MRDHKACLDP AKTQRLFFKW QQGHQLAIVN RFPGSKLTFQ RILSYDTVVQ TTGDSWHIPY DEDLVERCSC EVPTKWLCCLH QKT
<b>Source</b>	Baculovirus
<b>Target Names</b>	KIAA0391
<b>Protein Names</b>	Recommended name: Mitochondrial ribonuclease P protein 3 Short name= Mitochondrial RNase P protein 3
<b>Expression Region</b>	46-583
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.