



Recombinant Human Mitochondrial ribonuclease P protein 3 (KIAA0391)

Product Code	CSB-MP012187HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O15091
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	FSLKT MSPQNTKATN LIAKARYLRK DEGSNKQVYS VPHFFLAGAA KERSQMNSQT EDHALAPVRN TIQLPTQPLN SEEWDLKED LKENTGKTSF ESWIISQMAG CHSSIDVAKS LLAWVAAKNN GIVSYDLLVK YLYLCVFHMQ TSEVIDVFEI MKARYKTLEP RGYSLIRGL IHSDRWREAL LLEDIKKVI TPSKKNYNDQ IQGALLHQDV NTAWNLYQEL LGHDIVPMLE TLKAFFDFGK DIKDDNYSNK LLDILSYLRN NQLYPGESFA HSIKTWFESV PGKQWKGQFT TVRKSGQCSG CGKTIESIQL SPEEYECLKG KIMRDVIDGG DQYRKTPPQE LKRFFENFIKS RPPFDVVIDG LNVAKMFPKV RESQLLLNVV SQLAKRNLRL LVLGRKHMLR RSSQWSRDEM EEVQKQASCF FADDISEDDP FLLYATLHSG NHCRFITRDL MRDHKACLDP AKTQRLFFKW QQGHQLAIVN RFPGSKLTFQ RILSYDTVVQ TTGDSWHIPY DEDLVERCSC EVPTKWLCCLH QKT
Source	Mammalian cell
Target Names	KIAA0391
Protein Names	Recommended name: Mitochondrial ribonuclease P protein 3 Short name= Mitochondrial RNase P protein 3
Expression Region	46-583
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.