



# Recombinant Human WD repeat-containing protein 46 (WDR46)

<b>Product Code</b>	CSB-EP026034HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	O15213
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	METAPKPGKD VPPKKDKLQT KRKKPRRYWE EETVPTTAGA SPGPPRNKKN RELRPQRPN AYILKKSRI KKPQVPPKPR EWKNPESQRG LSGTQDPFPG PAPVPVEVVQ KFCRIDKSRK LPHSKAKTRS RLEVAEAE ETSIKAARSE LLLAEPEGFL EGEDGEDTAK ICQADIVEAV DIASAAKHFD LNLRQFGPYR LNYSRTGRHL AFGGRRGHVA ALDWVTKKLM CEINVMEAVR DIRFLHSEAL LAVAQNRWLH IYDNQGI ELH CIRRCDRVTR LEFLPFHLL ATASETGFLT YLDVSVGKIV AALNARAGRL DVMSQNPYNA VIHLGHSNGT VSLWSPAMKE PLAKILCHRG GVRVAVDST GTYMATSGLD HQLKIFDLRG TYQPLSTRTL PHGAGHLAFS QRGLLVAGMG DVVNIWAGQG KASPPSLEQP YLTHRLSGPV HGLQFCPFED VLGVGHTGGI TSMLVPGAGE PNFDGLESNP YRSRKQRQEW EVKALLEKVP AELICLDPRA LAEVDVISLE QGKKEQIERL GYDPQAKAPF QPKPKQKGRS STASLVKRKR KVMDEEHRDK VRQSLQQQHH KEAKAKPTGA RPSALDRFVR
<b>Source</b>	E.coli
<b>Target Names</b>	WDR46
<b>Protein Names</b>	Recommended name: WD repeat-containing protein 46 Alternative name(s): WD repeat-containing protein BING4
<b>Expression Region</b>	1-610
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.