



Recombinant Human WD repeat-containing protein 46 (WDR46)

Product Code	CSB-YP026034HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O15213
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	METAPKPGKD VPPKKDKLQT KRKKPRRYWE EETVPTTAGA SPGPPRNKKN RELRPQRPN AYILKKSRI KKPQVPPKPR EWKNPESQRG LSGTQDPFPG PAPVPVEVVQ KFCRIDKSRK LPHSKAKTRS RLEVAEAE EETSIIKAAARSE LLLAEPEGFL EGEDGEDTAK ICQADIVEAV DIASAAKHFD LNLRQFGPYR LNYSRTGRHL AFGGRRGHVA ALDWVTKKLM CEINVMEAVR DIRFLHSEAL LAVAQNRWLH IYDNQGIELH CIRRCDRVTR LEFLPFHLL ATASETGFLT YLDVSVGKIV AALNARAGRL DVMSQNPYNA VIHLGHSNGT VSLWSPAMKE PLAKILCHRG GVRVAVDST GTYMATSGLD HQLKIFDLRG TYQPLSTRTL PHGAGHLAFS QRGLLVAGMG DVVNIWAGQG KASPPSLEQP YLTHRLSGPV HGLQFCPFED VLGVGHTGGI TSMLVPGAGE PNFDGLESNP YRSRKQRQEW EVKALLEKVP AELICLDPRA LAEVDVISLE QGKKEQIERL GYDPQAKAPF QPKPKQKGRS STASLVKRKR KVMDEEHRDK VRQSLQQQHH KEAKAKPTGA RPSALDRFVR
Source	Yeast
Target Names	WDR46
Protein Names	Recommended name: WD repeat-containing protein 46 Alternative name(s): WD repeat-containing protein BING4
Expression Region	1-610
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.