



Recombinant *Saccharomyces cerevisiae* Maintenance of telomere capping protein 1 (MTC1)

Product Code	CSB-EP342449SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P47018
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MSENKNSEAE DVFEFLDPLP EAKNGGKMMVN TDVKGSQEGV KGGNSVAGK TGNDGKKGDD DIFEFLEELE KSNLSLTDKK GVEKKAPSES VNNKAQDEKV EESKENKNSE QDAHGKEKEP QQQEKEEEEE EEEEEEEEE ETPLHDPIAS ISNWWSSSGS AKVSSIWNKT AEQASQIKNR LAQEQLDLTS KINTSTITEI ARNLQKIVVG ETEEVLRHIL VHDLVNYPPL QYNIESKFDQ VLSSQVEGGI RIFVDEWGHP NNGGIPVEK KPSVADGELG NSKKKLQFNL FDGKVTDGEK LAFANLENAV KLFNTAHEEY QKQQKEADAT PDDDRSSISS NSNKISDLFI SILPIAIPQK QKDADGDFQV TDSNTPGNFN FTLVLKDITN DITITRSQG FPKVWVNWLE GSVEKTGSTA SEERNKSYDQ KKQKESEDED EDDEIIDPSE WVKEWIEDGL SLSFGVMAQN YVIDRMGL
Source	E.coli
Target Names	MTC1
Protein Names	Recommended name: Maintenance of telomere capping protein 1
Expression Region	1-478
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.