



Recombinant Arabidopsis thaliana Homeobox protein knotted-1-like 3 (KNAT3)

Product Code	CSB-EP342525DOA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P48000
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MAFHNNHLSQ DLSFNHFTDQ HQPPPPQPPP PPPQQQHFQ EAPPPNWLNT ALLRSSDNNN NFLNLHTATA NTTTASSSDS PSSAAAAAAA NQWLSRSSF LQRNNNNNAS IVGDGIDDVT GGADTMIQGE MKTGGENKN DGGGATAADG VVSWQNARHK AEILSHPLYE QLLSAHVACL RIATPVDQLP RIDAQLAQSQ HVVAKYSALG AAAQGLVGDD KELDQFMTHY VLLLCSFKEQ LQQHVRVHAM EAVMACWEIE QSLQSLTGVS PGEGMGATMS DDEDEQVESD ANMFDGGLDV LGFGPLIPTE SERSLMVRV QELKHEKQG YKEKIVDIRE EILRKRAGK LPGDTTSLVK AWWQSHSKWP YPTEEDKARL VQETGLQLKQ INNWFNQRK RNWHSNPSSS TVLKNKRKSN AGDNSGRERF A
Source	E.coli
Target Names	KNAT3
Protein Names	Recommended name: Homeobox protein knotted-1-like 3 Alternative name(s): Protein KNAT3
Expression Region	1-431
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.