



Recombinant Arabidopsis thaliana Aspartate carbamoyltransferase, chloroplastic (PYRB)

Product Code	CSB-MP342589DOA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P49077
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	CH AMQAGTRELK KFELSDVIEG KQFDREMLSA IFDVAREMEK IEKSSSQSEI LKGYLMTLF YEPSTRRLS FESAMKRLGG EVLTTENARE FSSAAKGETL EDTIRTVEGY SDIIVMRHFE SGAARKAAAT ANIPVINAGD GPGEHPTQAL LDVYTIQSEI GKLDGISVAL VGDLANGRTV RSLAYLLAKF KDVKIYFVSP EIVKMKDDIK DYLTSSGVEW EESSDLMEVA SKCDVVYQTR IQRERFGERL DLYEAARGKY IVDKDLLGVM QKKAIIMHPL PRLDEITADV DADPRAAYFR QAKNGLFIRM ALLKLLLVGW
Source	Mammalian cell
Target Names	PYRB
Protein Names	Recommended name: Aspartate carbamoyltransferase, chloroplastic EC= 2.1.3.2 Alternative name(s): Aspartate transcarbamylase Short name= ATCase
Expression Region	69-390
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Cytoplasmic domain
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.