



Recombinant Human Calpain-5 (CAPN5)

Product Code	CSB-BP004498HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O15484
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MFSCVKPYED QNYSALRRDC RRRKVLFPEDP LFPATDDSLY YKGTGPAVR WKRPKGICED PRLFVDGISS HDLHQGQVGN CWFVAACSSL ASRESLWQKV IPDWKEQEWD PEKPNAYAGI FHFHFWRFGE WVDVVIDDRL PTVNNQLIYC HSNSRNEFWC ALVEKAYAKL AGCYQALDGG NTADALVDFT GGVSEPIDLT EGDFANDETK RNQLFERMLK VHSRGGGLISA SIKAVTAADM EARLACGLVK GHAYAVTDVR KVRLGHGLLA FFKSEKLDMI RLRNPWGERE WNGPWSDTSE EWQKVKSER EKMGVTQDD GEFWMTFEDV CRYFTDIIKC RVINTSHLSI HKTWEEARLH GAWTLHEDPR QNRGGGCINH KDTFFQNPQY IFEVKKPEDE VLICIQRPK RSTRREGKGE NLAIGFDIYK VEENRQYRMH SLQHKAASSI YINSRSVFLR TDQPEGRYVI IPTTFEPGHT GEFLLRVFTD VPSNCRELRL DEPPHTCWSS LCGYPQLVTQ VHVLGAAGLK DSPTGANSYV IIKCEGDKVR SAVQKGTSTP EYNVKGIFYR KKLSQPITVQ VWNHRVLKDE FLGQVHLKAD PDNLQALHTL HLRDRNSRQP SNLPGTVAVH ILSSTSLMAV
Source	Baculovirus
Target Names	CAPN5
Protein Names	Recommended name: Calpain-5 EC= 3.4.22.- Alternative name(s): Calpain htra-3 New calpain 3 Short name= nCL-3
Expression Region	1-640
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Target Details	Calpains are calcium-dependent cysteine proteases involved in signal transduction in a variety of cellular processes. A functional calpain protein consists of an invariant small subunit and 1 of a family of large subunits. CAPN5 is one of the large subunits. Unlike some of the calpains, CAPN5 and CAPN6 lack a calmodulin-like domain IV. Because of the significant similarity to Caenorhabditis elegans sex determination gene tra-3, CAPN5 is also called as HTRA3.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.