



Recombinant Putative alpha-L-fucosidase (W03G11.3)

Product Code	CSB-MP342634CXY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P49713
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	DYTP DWESLDNRPL PSWYDDSKFG IFCHWGLYSV PAFRSEWMWW YWKGTQPKD VVNFVDKNIK PGTTYADFAK DFTAIFYNAN QFAETVKTSG ARYFVFTSKH HEGFTMWPSR TSWNWNMSMDI GPKRDIVGEL RDAFKKTDVH FGLYFSQFEW FHPMFLDDGK FNTTFYPEQV SYPQMIDIVT KYNPEVVWSD GEWDKSDDYW KAKEFLAWLY NSSPVKDQVV VNDRWGTGTM GKHGGFMTYS DHYDPGKLE KKWENCMTLD KHSWGNRRDM KASEVNTAYE IIEQLARTIA CNGNLLLNVG PNMHGQIPAI FEDRLEEIGR FVNITSEAIK GTRPWIHQND TSASNVWYTS KYSSGKKPLK NLYQNVYNFQ LEEHTIVYAW ILDTSHEQFE LKSVKTTKNT TATILGTDVV LTGFEESSDM IILSSKIDWK KLPRRDIIVL KIEKAASYLR NPLMSTNEHH VQ
Source	Mammalian cell
Target Names	W03G11.3
Protein Names	Recommended name: Putative alpha-L-fucosidase EC= 3.2.1.51 Alternative name(s): Alpha-L-fucoside fucohydrolase
Expression Region	17-482
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.