



Recombinant Zea mays Probable glutathione S-transferase BZ2 (BZ2)

Product Code	CSB-MP342684ZAX
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P50472
Product Type	Recombinant Protein
Immunogen Species	Zea mays (Maize)
Purity	≥85% (SDS-PAGE)
Sequence	MRVLGGEVSP FTARARLALD LRGVAYELLD EPLGPKKSDR LLAANPVYGK IPVLLLPDGR AICESAVIVQ YIEDVARESG GAEAGSLLLP DDPYERAMHR FWTAFIDDKF WPALDAVSLA PTPGARAQAA EDTRAALSLL EFAFKDRSNG RAFFSGGDAA PGLLDLALGC FLPALRACER LHGLSLIDAS ATPLLDGWSQ RFAAHPAAKR VLPDTEKVVQ FTRFLQVQAQ FRVHVS
Source	Mammalian cell
Target Names	BZ2
Protein Names	Recommended name: Probable glutathione S-transferase BZ2 EC= 2.5.1.18 Alternative name(s): Protein bronze-2
Expression Region	1-236
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.