



Recombinant Arabidopsis thaliana Glutathione S-transferase F2 (GSTF2)

Product Code	CSB-EP342801DOA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P46422
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	AGIKVFGHP ASIATRRVLI ALHEKNLDFE LVHVELKDGE HKKEPFLSRN PFGQVPAFED GDLKLFESRA ITQYIAHRYE NQGTNLLQTD SKNISQYAIM AIGMQVEDHQ FDPVASKLAF EQIFKSIYGL TTDEAVVAEE EAKLAKVLDV YEARLKEFKY LAGETFTLTD LHHIPAIQYL LGTPTKKLFT ERPRVNEWVA EITKRPASEK VQ
Source	E.coli
Target Names	GSTF2
Protein Names	Recommended name: Glutathione S-transferase F2 Short name= AtGSTF2 EC= 2.5.1.18 Alternative name(s): 24 kDa auxin-binding protein Short name= AtPM24 GST class-phi member 2
Expression Region	2-212
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.