



Recombinant Probable 26S protease regulatory subunit 6B (rpt-3)

Product Code	CSB-MP342806CXY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P46502
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MSAAALEEIG IAPLQENTSL RPPPLPATCV DVSSDDYIKL TTLERQLAHL QVMEDIKLE TRNLEKELLH AQEEVKRIQS VPLVIGQFLE AVDQNHAIVG STTGSNYYVR VLSILDRELL KPGCSVALHK YSNALVDVLP PEADSSIQML RPDEKPDISY GDIGGLDMQK QEVREAVELP LTHGELYQQI GIDPPRGVLM YGPPGCGKTM LAKAVAANTA ASFIRVVGSE FVQKYLGE GP RMVRDVFRLA KENSPSIIFI DEIDAIATKR FDAQTGADRE VQRILLELLN QMDGFDQSTN VKVIMATNRQ DTLDPALLRP GRLDRKIEFP LPDRRQKRLV FSTVCSRMLN SDDVDLEDWV ARPDKISGAD INSICQEAGM QAVRENRYVV LTKDLEKAYK NVVKKDTNDF EFYK
Source	Mammalian cell
Target Names	rpt-3
Protein Names	Recommended name: Probable 26S protease regulatory subunit 6B
Expression Region	1-414
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.