



Recombinant *Oryza sativa* subsp. japonica Homeobox protein knotted-1-like 6 (OSH1)

Product Code	CSB-YP342814OFG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P46609
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MEEISHHFGV VGASGVHGGH QHQHHHPWG SLSAIVAPP PPPQLQQQT QAGGMAHTPL TLNTAAAVG NPVLQLANGS LLDACGKAKE ASASASYAAD VEAIKAKIIS HPHYSSLAA YLDCQKVGAP PEVAARLTAV AQDLELRQRT ALGVLGAATE PELDQFMEAY HEMLVKYREE LTRPLQEAME FLRRVETQLN TLSISGRSLR NILSSGSSEE DQEGSGGETE LPEIDAHGVD QELKHLLKK YSGYLSSLKQ ELSKKKKKGK LPKDARQQLL NWWELHYKWP YPSESQKVAL AESTGLDLKQ INNWFINQRK RHWKPSDEMQ FVMMDGYHPT NAAFYMDGH FINDGGLYRL G
Source	Yeast
Target Names	OSH1
Protein Names	Recommended name: Homeobox protein knotted-1-like 6 Alternative name(s): Homeobox protein OSH1 Homeobox protein knotted-1-like 1 Short name= Oskn1
Expression Region	1-361
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.