



Recombinant *Saccharomyces cerevisiae* AP-3 complex subunit sigma (APS3)

Product Code	CSB-MP342854SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P47064
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MIHAVLIFNK KCQPRLVKFY TPVDLPKQKL LLEQVYELIS QRNSDFQSSF LVTPPSLLLS NENNNDEVNN EDIQIYKNY ATLYFTFIVD DQESELAILD LIQTFVESLD RCFTEVNELD LIFNWQTLES VLEEIVQGGM VIETNVNRIV ASVDELNCAA ESTDSKIGRL TSTGFGSALQ AFAQGGFAQW ATGQ
Source	Mammalian cell
Target Names	APS3
Protein Names	Recommended name: AP-3 complex subunit sigma Alternative name(s): AP-3 complex sigma3A subunit Adapter-related protein complex 3 sigma subunit Clathrin-associated/assembly/adaptor protein, small 3 Sigma-adaptin 3A Sigma3-adaptin
Expression Region	1-194
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.