



Recombinant *Saccharomyces cerevisiae* Uncharacterized protein YJR111C (YJR111C)

Product Code	CSB-EP342863SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P47148
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MNQILNAQRL IQLSQFHPKL KNIWYLVAAA TFSVCNEPQE IPKLYHYAML LSNDNAHMYR FTLASQTIDL LRSELP MRKT LINENYQQPT FFQKQLTAKF REVILKTGPL AGLPRAINGL TMLKETT PDI LVP HLDPIDP WEAAMGNSSP LSETSMRRKH DKTIQERDHT IQNGLRHWNS IYNKVSTRVV>NNLNSSYPDL WYYTLVHVYG PLFAFDEILS AQETSLVIIA SLVPQDVNPQ LRGHLKGALN IGCDKETVEA VRGLAILISQ WCGVSWKSGV VKL
Source	E.coli
Target Names	YJR111C
Protein Names	Recommended name: Uncharacterized protein YJR111C
Expression Region	1-283
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.