



Recombinant Putative serine/threonine-protein phosphatase C23G10.1 (C23G10.1)

Product Code	CSB-MP342945CXY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P48459
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MSKKPRVPEV AEMSAVSDND VIKEKGARDK MFQKQSDESN KMFAEHFIKT LLACKGMTKI RTMDIFRLIH ICKKIFTVQK SMVEIDGPVR ICGDLHGQYP DLIRLFAQVL GGFPPDSNYL FLGDYVDRGS FNLEVILLCL AYKARYPNNF MMLRGNHEVI HINEKYGFKD EVFNRRKGEYH DELYPEFNEM MDMMPLVALV GGRILCMHGG LSQHIKSLDD LRNLRRPFHS EDECLENDIM WSDPAKVSGW TANPRGASVQ FGENEVKEMC KLLDIDLIVR GHQVVQDGYE FFAGKKLVTV FSAPHYMQSF TNSAAVCKVS AGLEVSFEVL KPEDIRVEEI KCSAESSCAS DMQQ
Source	Mammalian cell
Target Names	C23G10.1
Protein Names	Recommended name: Putative serine/threonine-protein phosphatase C23G10.1 EC= 3.1.3.16
Expression Region	1-354
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.