



Recombinant *Saccharomyces cerevisiae* Ubiquitin domain-containing protein DSK2 (DSK2)

Product Code	CSB-EP342949SVG-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P48510
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MSLNIIHKSG QDKWEVNVAP ESTVLQFKEA INKANGIPVA NQRLIYSGKI LKDDQTVESY HIQDGHSVHL VKSQPKPQTA SAAGANNATA TGAAAGTGAT PNMSSGQSAG FNPLADL TSA RYAGYLNMP S ADMFGPDGGA LNNSDNNQDE LLRMMENPIF QSQMNEMLSN PQMLDFMIQS NPQLQAMGPQ ARQMLQSPMF RQMLTNPDMI RQSMQFARMM DPNAGMGSAG GAASAFPAPG GDAPEEGSNT NTTSSSNTGN NAGTNAGTNA GANTAANPFA SLLNPALNPF ANAGNAASTG MPAFDPALLA SMFQPPVQAS QAEDTRPPEE RYEHQLRQLN DMGFFDFDRN VAALRRSGGS VQGALDSSLN GDV
Source	E.coli
Target Names	DSK2
Protein Names	Recommended name: Ubiquitin domain-containing protein DSK2
Expression Region	1-373
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.