



Recombinant Probable protein phosphatase 2C T23F11.1 (T23F11.1)

Product Code	CSB-BP343023CXY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P49596
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	MGQTLSEPVT KKEASCAN E NYLVGSSCMQ GWRVDMEDAH THLLSLPDDP KCAFFAVYDG HGGSKVSQYS GINLHKVVA QKEFSEGNMK EAIEKGFLEL DQQMRVDEET KDDVSGTTAV VVLIKEGDVY CGNAGDSRAV SSVVGEARPL SFDHKPSHET EARRIIAAGG WVEFN RVNGN LALSRLGDF AFKNCDTKPA EEQIVTAFPD VITDKLTPDH EFIVLACDGI WDVMTNQEVV DFVREKLAEK RDPQSICEEL LTRCLAPDCQ MGGLGCDNMT VVLVGLLHGQ SPDTLFTKCA RPAVFTGVNN EDGDQNIQQ DISRVINFD GEQRVNAANQ EEEEDDNEPA PANFQV
Source	Baculovirus
Target Names	ppm-2
Protein Names	Recommended name: Probable protein phosphatase 2C T23F11.1 Short name= PP2C EC= 3.1.3.16
Expression Region	1-356
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.