



Recombinant *Carassius auratus* Transcriptional regulator Myc (myc)

Product Code	CSB-YP343029DQG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P49709
Product Type	Recombinant Protein
Immunogen Species	<i>Carassius auratus</i> (Goldfish)
Purity	>85% (SDS-PAGE)
Sequence	MPVSASLAYK NYDYDYDSIQ PYFYFDNDDE DFYHHQQGQP QPPAPSEDIW KKFELLPTPP LSPSRRQSLs TAEQLEMVSE FLGDDVvNQS FICDADYSQS FIKSIIIQDC MWSGFSAAAK LEKAVSERLA SLHAARKELI SDSSSNRLSA SYLQDLSTSA SECIDPSVVF PYPLTESSKS NKVAPSQPML VLDTPPNSSS SSGSDSEDEE EEEEEEEEEEE EEEEEEEEEEE EIDVVTVEKR QKRNEADVSD SRYPSPVLK RCHVSTHQHN YAAHPSTRHD QPAVKRLRLE TSSSNRHGKQ RKCTSPRTSD SEDNDKRRTH NVLERQRRNE LKLSFFALRD EIPEVANNEK AAKVVILKKA TECIHSMQLD EQRLLSIKEQ LRRKSEQLKH RLQQLRSSH
Source	Yeast
Target Names	myc
Protein Names	Recommended name: Transcriptional regulator Myc Short name= c-Myc
Expression Region	1-399
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.