



# Recombinant Agropine synthesis reductase (mas1)

<b>Product Code</b>	CSB-EP343063AEW-B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P50201
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Agrobacterium rhizogenes
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MPISLEYLGP AQAPASVGDY YFLRHGKTDL NEKEVVFVQGE KHWGVQAGAGT NIGLNETGKR QAVLAGNVLR KLPISVVCS PLLRAIQTAV IANIGCLCFE IDDDLKERDF GKHEGGYGPL KMFEENYPDT ENTELFMRV AKALAAHAE NTLFVSHGGV LRVVAALLGV ELTKEHTDNG RVLHFRRGCP QWTVEIHQSP VILISGPSRG IGKAI AENLI AHGYRMSLGA RRVKDLESAF GPQSEWLHYA RFDAEDNDTM AAWVTA AVEK FGRIDGLINN AGCGEPV NLE KDIDYKQFHR QWHINCVAPL RMTCLCLPYL CETGSGRVVN INSMGQ RVL NPFVGYNMTK HALGGLTKTT QHVGWDRGVR AIDICPGFVA TEMSAWTDLI SSNDMIQPND IANLVREAIE RPNRAYVPRS EVMCMKESIH
<b>Source</b>	E.coli
<b>Target Names</b>	mas1
<b>Protein Names</b>	Recommended name: Agropine synthesis reductase EC= 1.-.-.
<b>Expression Region</b>	1-430
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.