



Recombinant Mouse Serine/threonine-protein kinase Chk1 (Chek1)

Product Code	CSB-MP005338MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O35280
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MAVPFVEDWD LVQTLGEGAY GEVQLAVNRI TEEAVAVKIV DMKRAIDCPE NIKKEICINK MLSHENVVKF YGHRREGHIQ YLFLEYCSGG ELFDRIEPTD GMPEQDAQRF FHQLMAGVVY LHGIGITHRD IKPENLLLDE RDNLKISDFG LATVFRHNNR ERLLNKMCMT LPYVAPPELLK RKEFHAEPVD VWSCGIVLTA MLAGELPWDQ PSDSCQEYSD WKEKKTYLNP WKKIDSAPLA LLHKILVETP SARITIPDIK KDRWYNKPLN RGAKRPRATS GGMSESSSGF SKHIHSNLDF SPVNGSSEE TVKFSSSQPE PRTGLSLWDT GPSNVDKLVQ GISFSQPTCP EHMLVNSQLL GTPGSSQNPW QRLVKRMTRF FTKLDADKSY QCLKETFEKL GYQWKKSCMN QTVSTTDRR NNKLIFKINL VEMDEKILVD FRLSKGDGLE FKRHFLKIKG KLSDVVSSQK VWFPVT
Source	Mammalian cell
Target Names	Chk1
Protein Names	Recommended name: Serine/threonine-protein kinase Chk1 EC= 2.7.11.1 Alternative name(s): CHK1 checkpoint homolog Checkpoint kinase-1
Expression Region	1-476
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.