



Recombinant *Oryza sativa* subsp. japonica 26S protease regulatory subunit 4 homolog (TBP2)

Product Code	CSB-MP343201OFG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P46466
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MGQGTPGGMG KQGGAPGDRK PGGDGDKKDR KFEPAAAPSR VGRKQRKQKG PEAAARLPNV APLSKCRLRL LKLERVKDYL LMEEEFVAAQ ERLRPTEDKT EEDRSKVDDL RGTPMSVGS L EEIIDESHAI VSSSVGPEYY VGILSFVDKD QLEPGCSILM HNKVLSVVG I LQDEVDPMS VMKVEKAPLE SYADIGGLDA QIQEIKEAVE LPLTHPELYE DIGIRPPKGV ILYGEPGTGK TLLAKAVANS TSATFLRVVG SELIQKYLGD GPKLVRELF R VADELSPSIV FIDEIDAVGT KRYDAHSGGE REIQRMTLEL LNQLDGFDSR GDVKVILATN RIESLDPALL RPGRIDRKIE FPLPDIKTRR RIFQIHTSKM TLADDVNLEE FVMTKDEFSG ADIKAICTEA GLLALRERRM KVTHADFKKA KEKVMFKKKE GVPEGLYM
Source	Mammalian cell
Target Names	TBP2
Protein Names	Recommended name: 26S protease regulatory subunit 4 homolog Alternative name(s): Tat-binding protein homolog 2
Expression Region	1-448
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.