



# Recombinant SAM50-like protein gop-3 (gop-3)

<b>Product Code</b>	CSB-EP343211CXY
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P46576
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Caenorhabditis elegans
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSEKTFHKAQ TIRAKASGVP SIVEAVQFHG VRITKNDALV KEVSELYRSK NLDELVHNSH LAARHLQEVG LMDNAVALID TSPSSNEGYV VNFLVREPKS FTAGVKAGVS TNGDADVSLN AGKQSVGGRG EAINTQYTYT VKGDHCFNIS AIKPFLGWQK YSNVSATLYR SLAHMPWNQS DVDENAVALA YNGQLWNQKL LHQVKLNAIW RTLRA TRDAA FSVREQAGHT LKFSLENAVA VDTRDRPILA SRGILARFAQ EYAGVFGDAS FVKNTLDLQA AAPLPLGFIL AASFQAKHLK GLGDREVVHL DRCYLGGQQD VRGFG LNTIG VKADNSCLGG GASLAGVVHL YRPLIPP NML FAHAFLASGS VASVH SKNLV QQLQDTQRVS AGFGLAFVFK SIFRLELNYT YPLKYVLGDS LLGGFHIGAG VNFL
<b>Source</b>	E.coli
<b>Target Names</b>	gop-3
<b>Protein Names</b>	Recommended name: SAM50-like protein gop-3 Alternative name(s): Gro-1 operon protein 3
<b>Expression Region</b>	1-434
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.