



Recombinant Plasmodium falciparum Plasmepsin-2

Product Code	CSB-EP343237PLO-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P46925
Product Type	Recombinant Protein
Immunogen Species	Plasmodium falciparum
Purity	>85% (SDS-PAGE)
Sequence	SSNDNI ELVDFQNIMF YGDAEVGDNQ QPFTFILD TG SANLWVPSVK CTTAGCLTKH LYDSSKSR TY EKDGTKVEMN YVSGTVSGFF SKDLVTVGNL SLPYKFIEVI DTNGFEPTYT ASTFDGILGL GWKDLSIGSV DPIVVELKNQ NKIENALFTF YLPVHDKHTG FLTIGGIEER FYEGPLTYEK LNHDLYWQIT LDAHVGNIML EKANCIVDSG TSAITVPTDF LNKMLQNL DV IKVPFLPFYV TLCNNSKLPT FEFTSENGKY TLEPEYYLQH IEDVGPGLCM LNIIGLDFPV PTFILGDPFM RKYFTVFDYD NHSVGIALAK KNL
Source	E.coli
Protein Names	Recommended name: Plasmepsin-2 EC= 3.4.23.39 Alternative name(s): Aspartic hemoglobinase II PFAPD
Expression Region	125-453
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.